



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Termopar de unin expuesta IEC de un solo
disparo de PTFE con enchufe en miniatura - Tipo
T - Longitud de cable de 300 mm**

Image
Coming
Soon

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0439-001
General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
Sensor Type	Exposed Junction
Thermocouple Type	T
Cable Length	300mm
Cable Type	PTFE
Core / Strands	1/0.376mm
Termination Type	Terminated with a miniature flat pin plug (Plug max temp 220°C)
Max. Temperature	250°C
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Colour	Brown
Series	Exposed Junctions