



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Termopar de unin expuesta aislado con fibra de vidrio ANSI con enchufe en miniatura - Tipo K - 300 mm de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-0457-001
General Description	Fast response, welded exposed junction
Sensor Type	Good temperature resistance (not suitable in fluids)
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	300mm
Cable Type	glassfibre insulated flat pair cable
Core / Strands	1/0.3mm conductors
Cable Termination	Miniature Plug
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	ANSI MC96.1
Connector Type	Plug
Connector Size	Miniature
Colour	Yellow
Series	Exposed Junctions