


**Termopar de unin expuesta aislado con PFA ANSI
con enchufe en miniatura - Tipo K - Tamao de
hebra de 1 / 0.3 mm - Longitud de 300 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-0454-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	300mm
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI MC96.1
Series	Exposed Junctions