




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de unin expuesta aislado ANSI PFA con colas desnudas - Tipo K - Tamao de hebra de 1 / 0.3 mm - Longitud de 300 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

<b>Product Code</b>	XE-0434-001
<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	K
<b>Thermocouple Type</b>	Exposed Junction
<b>Cable Length</b>	300mm
<b>Cable Type</b>	Teflon® PFA insulated twin-twisted lead
<b>Core / Strands</b>	1/0.3mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI MC96.1
<b>Series</b>	Exposed Junctions