


## Conector de termopar de cable rpido en miniatura AM-N-FQ Tipo N ANSI



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de cable rpido en miniatura para toma de termopar AM-N-FQ Tipo N ANSI**

Conexin fcil con el pinchn, solo hay que empujar el cable y apretar el tornillo. Los contactos estn polarizados para evitar terminaciones incorrectas.

**Especificaciones****Specifications**

|                            |                                                            |
|----------------------------|------------------------------------------------------------|
| <b>Product Code</b>        | XE-1282-001                                                |
| <b>General Description</b> | Easy jab-in connection, just push in wire & tighten screw. |
| <b>Thermocouple Type</b>   | N ANSI                                                     |
| <b>Connector Type</b>      | Quick Wire Socket                                          |
| <b>Connector Size</b>      | Miniature                                                  |
| <b>Colour</b>              | Orange                                                     |
| <b>Max. Temperature</b>    | 220°C                                                      |
| <b>+Positive Contact</b>   | Nickel / Chromium / Silicon                                |
| <b>-Negative Contact</b>   | Nickel Silicon                                             |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                                     |