


Conector de termopar de alambre rpido en miniatura AM-R/S-FQ Tipo R/S ANSI



Miniature in-line Socket - Quick Wire
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de cable rpido en miniatura Toma de termopar AM-R/S-FQ Tipo R/S ANSI

Conexin fcil con el pinchn, solo hay que empujar el cable y apretar el tornillo. Los contactos estn polarizados para evitar terminaciones incorrectas.

Especificaciones**Specifications**

Product Code	XE-1322-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	R/S ANSI
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS