

Sonda de carne in almblica RFX



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockists of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Realice un seguimiento de las temperaturas centrales de la carne desde su dispositivo

Utilice hasta 50 sondas con un GATEWAY

Conexión inalámbrica fiable

50 horas de duración de la batería con una carga de 10 minutos

Totalmente resistente al agua hasta IP69K

Certificado de calibración trazable gratuito

Debe estar emparejado con
0

[RFX Gateway](#)
[RFX Starter Kit](#)

Características técnicas

Specifications

Product Code	XE-4201-001
Model	RFX Wireless Meat Probe
Type	Thermistor, 4 Sensors
Accuracy	+/-0.5°C (-10 to 100°C)
Battery	1.5 V AAA with charging block - 175+ full charges per AAA battery. 0-80% in 10 minutes; 0-100% in 20 minutes
warranty	6 Months
Output	Wireless communication: 433Mhz radio frequency
Temperature Range	Internal: -18 to 100°C / External: -18 to 538°C
Standards Met	IP69K
Measuring range	Battery powered RFX GATEWAY: 650m (2,132 ft.), AC Powered RFX GATEWAY: 385m (1,263 ft.) (RFX MEAT to RFX GATEWAY)