

Sondas de temperatura de termopar atornilladas ANSI - Tipo K - 1/4" BSPP - 100mm de longitud x 6mm de diametro



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sensores de temperatura tipo K atornillados con un rango de medicin de -75 a +250 C y cable conductor de par plano PFA de 1 metro.

Para montaje directo en recipientes, tuberías, tanques o conductos para la medicin de temperatura de líquidos o gases.

Aplicaciones tpicas

Control de procesos industriales

Tuberías o conductos

Recipientes y tanques

Caracterstias tcnicas

Specifications

Product Code	XE-7369-001
General Description	Screw in Type K thermocouple
Thermocouple Type	K ANSI
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	100mm
Sheath Diameter (mm)	6mm
Process Connection options	1/4"BSPP
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI
Number of Wires	2