


Termopar de unin expuesta IEC con doble enlace aislado en PTFE con tapn miniatura - Tipo K - 3 M de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar tipo K. Respuesta rpida, unin soldada expuesta, disponible en la especificacin de cable Clase 1 conforme a la BS EN 60584-3:2008 con terminacin de enchufe miniatura.

Ideal para aplicaciones de prueba y desarrollo.

Especificaciones

Specifications

Product Code	XE-3507-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	3M
Cable Type	Teflon® PTFE insulated twin-twisted lead
Core / Strands	1/0.2mm
Termination Type	Miniature flat pin plug colour coded in accordance IEC 584
Max. Temperature	260°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Series	Exposed Junctions