


# Conector de termopar para montaje de PCB AM-T-PCB Tipo T ANSI - Pines chapados en estaño



**Miniature PCB Socket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar de montaje de PCB Socket AM-T-PCB-TIN Tipo T ANSI

Especialmente diseñados para PCB, pines de estao macizos para montaje directo en placa y provisin para sensor CJC

### Especificaciones

## Specifications

<b>Product Code</b>	XE-1230-001
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	T ANSI
<b>Connector Size</b>	Miniature
<b>Connector Type</b>	PCB socket
<b>Colour</b>	Blue
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Tin Plated
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS