


Conector de termopar para montaje en miniatura de panel de cable rpido tipo J ANSI



Miniature Socket -

Quick Wire with stainless steel mounting bracket

(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Montaje de panel para conector de termopar de alambre rpido en miniatura con soporte de acero inoxidable
AM-J-SSPFQ Tipo J ANSI**

Portacasillas de tamao miniatura suministradas con soportes de acero inoxidable, permiten un ajuste rpido, fiable y resistente en el panel

Especificaciones**Specifications**

Product Code	XE-9066-001
General Description	Contacts are polarized to ensure correct connection
Thermocouple Type	J ANSI
Connector Type	Quick Wire Socket with SS Bracket
Connector Size	Miniature
Colour	Black
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS