

Extensores de sondas EL-LASCAR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

EXTENSORES DE SONDAS EL

Jack de 3,5 mm a zcalo hembra de 3,5 mm que puede usarse para extender la sonda suministrada con todos los productos Thermistor (TP y TP+) de Lascar.

Especificaciones**Specifications**

General Description	Probe extenders for Lascar data loggers
Type	3.5mm jack to female 3.5mm socket
Cable Length	1.5M, 3M, 5M, 10M