


# Conector de termopar de montaje en miniatura de panel de alambre rpido para soporte de acero inoxidable AM-E-SSPFQ Tipo E ANSI



**Miniature Socket -**

**Quick Wire with stainless steel mounting bracket**

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Montura de panel para conector de termopar de alambre rpido en miniatura con soporte de acero inoxidable  
AM-E-SSPFQ Tipo E ANSI**

Portacasillas de tamao miniatura suministradas con soportes de acero inoxidable, permiten un ajuste rpido, fiable y resistente en el panel

**Especificaciones****Specifications**

<b>Product Code</b>	XE-9036-001
<b>General Description</b>	Contacts are polarized to ensure correct connection
<b>Thermocouple Type</b>	E IEC
<b>Connector Type</b>	Quick Wire Socket with SS Bracket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Violet
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS