




Conector de termopar para montaje de PCB Socket AM-K-PCB Tipo K ANSI



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar para montaje de PCB Socket AM-K-PCB Tipo K ANSI

Diseado especialmente para PCB, pines slidos para montaje directo en placa y provisin para sensor CJC

Especificaciones**Specifications**

Product Code	XE-1139-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	K ANSI
Connector Type	PCB socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector