


Conector de termopar miniatura JM-T-M Tipo T JIS



Miniature in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar miniatura JM-T-M Tipo T JIS

Con pines planos adaptados al tipo termopar y ideales para uso general, todos los contactos estn polarizados para asegurar una conexin correcta. Estos conectores aceptan cables de termopar de hasta 4 mm de dimetro.

Especificaciones**Specifications**

Product Code	XE-9250-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	T JIS
Connector Type	Plug
Connector Size	Miniature
Colour	Brown
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector