


## Conector de termopar de alambre rpido en miniatura JM-T-FQ Tipo T JIS



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de cable rpido miniatura para toma de termopar JM-T-FQ Tipo T JIS**

Conexin fcil con el pinchn, solo hay que empujar el cable y apretar el tornillo. Los contactos estn polarizados para evitar terminaciones incorrectas.

**Especificaciones****Specifications**

<b>Product Code</b>	XE-9263-001
<b>General Description</b>	Easy jab-in connection, just push in wire & tighten screw.
<b>Thermocouple Type</b>	T JIS
<b>Connector Type</b>	Quick Wire Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Brown
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS