


## Conector de termopar de vaso rectangular en miniatura JM-J-FF Tipo J



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar rectangular miniatura Socket de fascia JM-J-FF Tipo J**

Pinza de nylon para fijar rpido y sencillo al panel. Todos los zcalos de montaje de la fascia estn polarizados mediante pines de tamao desigual

**Especificaciones****Specifications**

<b>Product Code</b>	XE-9177-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	J JIS
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Iron
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS