


Conector termopar estndar de montaje en panel con soporte de acero inoxidable JS-T-SSPF Tipo T JIS



Standard in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector termopar estndar de montaje en panel con soporte de acero inoxidable JS-T-SSPF Tipo T JIS

Pensado para su uso en procesos ms industriales, pero puede emplearse en todas las aplicaciones segn se desee

Especificaciones**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9378-001 |
| General Description | Intended for use in more industrial processes but can be used in all applications as desired |
| Thermocouple Type | T JIS |
| Connector Type | Socket with SS Bracket |
| Connector Size | Standard |
| Colour | Brown |
| Max. Temperature | 220°C |
| +Positive Contact | Copper |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |