


Conector estndar de termopar de montaje rpido en panel con soporte de acero inoxidable JS-B-SSPFQ Tipo B JIS



Standard Socket - Quick Wire with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector estndar de termopar de montaje rpido en panel con soporte de acero inoxidable JS-B-SSPFQ Tipo B JIS

Pensado para su uso en procesos ms industriales, pero puede emplearse en todas las aplicaciones segn se desee.
Conexin sencilla de 'jab-in' para terminacin rpida, solo hay que empujar el cable y apretar el tornillo.

Especificaciones**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-9339-001 |
| General Description | Easy 'jab-in' connection for rapid termination, just push wire in and tighten screw |
| Thermocouple Type | B JIS |
| Connector Type | Quick Wire Socket with SS Bracket |
| Connector Size | Standard |
| Colour | Grey |
| Dimensions | See Schematic Drawing (Product Literature Tab) |
| Max. Temperature | 220°C |
| +Positive Contact | Copper |
| -Negative Contact | Copper |
| Standards Met | BS EN 50212:1997. RoHS |